

# AUTOMATED FILTRATION AND WASHING PLATFORM

## Why study Filtration and Washing

- Isolation and drying can account for more than 50% of process duration
- Identify crystallization and isolation issues at an early stage - helps design of better processes
- Optimise wash protocol (quantities and time) to reduce cycle time and improve product quality
- Gain knowledge of filtration and washing time - helps in plant occupancy planning

## Automated Filtration and Washing Platform



- Bespoke automated self-contained platforms installed in pharmaceutical industry
- Utilises imaging/vision technology for intelligent feedback and real time control
- Data systematically gathered and stored in secured database, with automated report generation

## What does it do?

### Filtration & Washing of...

- Up to 50ml slurry
- Up to 6 washes, with different wash composition
- Multiple samples

### Variable Parameters:

- Pressure (0.2-1.8 Bar)
- Temperature of experiment (10-50°C)
- End-point selection - dryland or breakthrough
- Up to 15 solvent choices
- Pause for cake sampling
- Contact and "additional" time


### Data collections:

- Setting test
- Online viscosity
- Filtrate mass with time
- Cake volume and depth
- Automated sampling and dilution of filtrates for analysis
- Segregation and weighing of all filtrate and retains filter cake

### Health and Safety

- Fully contained platform
- Automatic shutdown in events of ventilation failure
- CE and UL marked

## Images from the Platform




- Images of the filter cake and filtrate are captured
- Help understanding poor wash efficiency
- Early warning of tendency of cake cracking or pulling from the sides

## Automated Reporting System


- Searchable database with 3 types of reports

Sample ID	Report	Start Date	End Date	Pressure (Bar)	Temp (°C)	Filter Unit	Comments
1	AFnW	20-01-2008	20-01-2008	1.8 bar	20	100000	100000
2	AFnW	20-01-2008	20-01-2008	1.8 bar	20	100000	100000
3	AFnW	20-01-2008	20-01-2008	1.8 bar	20	100000	100000
4	AFnW	20-01-2008	20-01-2008	1.8 bar	20	100000	100000
5	AFnW	20-01-2008	20-01-2008	1.8 bar	20	100000	100000

- Data manipulation is based on Darcy's law; cake resistance and cake permeability are calculated automatically; PDF file generated; option to combine data to provide cake compressibility coefficient report



- Advanced report gives all additional data collected including calculations



## Where AFnW adds value

- Filtration under controlled identical conditions to
  - Compare batches from different crystallization conditions
  - Compare batches prepared using different input source
  - Determine the effect of scale-up on product formed. (Lab versus scale-up versus pilot plant batches)
  - Benchmark plant equipment
- Predict filtration on scale up
- Determine optimum filtration and washing conditions
- Troubleshooting
  - Sub-sampling from plant to compare actual versus proposed conditions

## Other Potential Uses

- Isolation and efficient washing of catalyst
- Removal of inorganics
- Characterisation of the removal of silica, alumina or charcoal from organic synthetic mixtures

## Example # 1: Scale Up Calculations

- Cake resistance data can be used to select optimum process conditions and calculate filtration time at scale





- Plant scale-up predicted filtration time at 0.5bar (20kg batch, 0.2m<sup>2</sup>):
  - 2 hours for non seeded crystallisation
  - 20min for seeded crystallisation
  - The range of pressure was investigated (0.5-1.5bar) - low compressibility index
- 1.5bar selected,
  - 15min expected filtration time on scale,
  - 15-30min filtration time on plant observed

## Example # 2: Wash Efficiency Curves

- Wash efficiency curve gives an indication of the number of washes required to control key impurities

## Wash Volumes Used for Wash Efficiency Curve

Original Process Volume	1 <sup>st</sup> Process Volume	2 <sup>nd</sup> Process Volume	3 <sup>rd</sup> Process Volume	4 <sup>th</sup> Process Volume	5 <sup>th</sup> Process Volume	6 <sup>th</sup> Process Volume
Wash Volume on AFnW	1 Cake Volume	2 Cake Volumes	3 Cake Volumes	4 Cake Volumes	5 Cake Volumes	6 Cake Volumes
Mother Liquor	Wash 1	Wash 2	Wash 3	Wash 4	Wash 5	Wash 6



- Results show good removal of impurities during isolation process with 2 process volumes
- Wash volumes can be reduced from 3 to 2 process volumes