

Specifications

Description **Stirrer Assembly for pressure and vacuum applications**

Part number	B19	20080
	B24	20160
	B29	20170

Conformity		CE
Commodity Code		3917290099

Materials of construction

Body		PEEK
Seal		Spring energised seal
Shaft diameter	mm	10mm
Max operating temperature	°C	100
Max operating pressure	bar	10

Shaft material options		Precision Glass 316 Stainless Steel
------------------------	--	--

Environment

Permissible ambient temperature	°C	5 - 40
---------------------------------	----	--------

Rodaviss joints and connectors

<http://www.duran-group.com/en/products-solutions/glassware-for-manipulators/products/joints- connectors/rodaviss.html>