

Specifications

Description **Stirrer Assembly for pressure and vacuum applications**

Conformity CE
Commodity Code 3917290099

Part number	B14	51890
	B19	83260
	B24	84330
	B29	84360

Materials of construction

Body PEEK

Seal Spring energised Teflon seal

Spring Material 316 Stainless Steel

Shaft diameter	B14	mm	6
	B19	mm	10
	B24	mm	10
	B29	mm	10

Max operating temperature °C 80

Max operating pressure bar 10

Shaft material options Precision Glass
316 Stainless Steel

Environment

Permissible ambient temper: °C 5 - 40

Rodaviss joints and connectors

<http://www.duran-group.com/en/products-solutions/glassware-for-manipulators/products/joints-connectors/rodaviss.html>